

Análisis transcriptómicos de la expresión génica

Máster Universitario en Bioinformática

Sesión 7



Universidad
Internacional
de Valencia

Dra. Paula Soler Vila
paula.solerv@professor.universidadviu.com

De:
 Planeta Formación y Universidades




Bloque III: Análisis de datos de NGS




Objetivos de la sesión

- 1 Comprender las principales características de los **alineadores** y ejecutar el proceso de alineamiento.
- 2 Analizar y comprender la información contenida en los archivos de alineamiento en formato **SAM/BAM**.
- 3 Analizar y comprender la información contenida en los archivos de anotaciones en formato **GFF/GTF**.
- 4 Realizar la **cuantificación** de los alineamientos utilizando un archivo de anotaciones como referencia.

Flujo de trabajo del análisis de datos de RNA-seq (NGS)



HISAT2 (Hierarchical Indexing for Spliced Alignments of Transcripts)



HISAT2

graph-based alignment of next generation sequencing reads to a population of genomes

HISAT2 is a fast and sensitive alignment program for mapping next-generation sequencing reads (both DNA and RNA) to a population of human genomes as well as to a single reference genome. Based on an extension of BWT for graphs (Sirén et al. 2014), we designed and implemented a graph FM index (GFM), an original approach and its first implementation. In addition to using one global GFM index that represents a population of human genomes, **HISAT2** uses a large set of small GFM indexes that collectively cover the whole genome. These small indexes (called local indexes), combined with several alignment strategies, enable rapid and accurate alignment of sequencing reads. This new indexing scheme is called a Hierarchical Graph FM index (HGFM).

The HISAT-3N paper published at *Genome Research*. 7/1/2021

HISAT-3N beta release 12/14/2020

HISAT-3N is a software system for analyzing nucleotide conversion sequencing reads. See the [HISAT-3N](#) for more details.

[Main](#)

[About](#)

[Manual](#)

[HISAT-3N](#)


[Download](#)

[HowTo](#)

[Links](#)

<http://daehwankimlab.github.io/hisat2/>

HISAT2 -> Opciones mínimas a considerar



HISAT2

graph-based alignment of next generation sequencing reads to a population of genomes

- **-x** : **prefijo** del índice del genoma de referencia [**genome**]
- **-U** : lista de lecturas para ser alineadas [**trimmed**]
- **-S** : archivo de salida en formato SAM
- **-k** : define el número máximo de alineamientos por lectura.

PRACTIQUEMOS

```
(05MBIF) Data]$ tree -L 4
├── Annotation
├── Processed
│   ├── 01.Quality_control
│   │   ├── SRR1552444_fastqc.html
│   │   └── SRR1552444_fastqc.zip
│   ├── 02.Trimming
│   │   ├── SRR1552444.fastq.gz_trimming_report.txt
│   │   └── SRR1552444_trimmed.fq.gz
│   └── 03.Alignment
│       └── SRR1662444_hisat2.sam
├── Raw
│   └── SRR1552444.fastq.gz
├── Reference_genome
│   └── mm10
│       ├── genome.1.ht2
│       ├── genome.2.ht2
│       ├── genome.3.ht2
│       ├── genome.4.ht2
│       ├── genome.5.ht2
│       ├── genome.6.ht2
│       ├── genome.7.ht2
│       └── genome.8.ht2
│       ├── make_mm10.sh
│       └── mm10_genome.tar.gz
```



HISAT2 -> Realizando el alineamiento




```
(05MBIF) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr 03.Alignment]$
```


```
hisat2 -k1 -U ../02.Trimming/SRR1552444_trimmed.fq.gz -x ../Reference_genome/mm10/  
genome -S SRR1552444_hisat2.sam
```

- **-x** : **prefijo** del índice del genoma de referencia [*genome*]
- **-U** : lista de lecturas para ser alineadas [*trimmed*]
- **-S** : archivo de salida en formato SAM
- **-k** : define el número máximo de alineamientos por lectura.

Flujo de trabajo del análisis de datos de RNA-seq (NGS)



Flujo de trabajo del análisis de datos de RNA-seq (NGS)



HISAT2 -> Realizando el alineamiento -> Resultado

27906762 reads; of these:

27906762 (100.00%) were unpaired; of these:

817651 (2.93%) aligned 0 times

27089111 (97.07%) aligned exactly 1 time

0 (0.00%) aligned >1 times

97.07% overall alignment rate

El porcentaje debería ser mayor del 80%

-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler users 8.3G Jul 27 11:02 SRR1552444_hisat2.sam

Archivo de alineamientos SAM

```
(05MBIF) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr 03.Alineamiento]$ head SRR1552444_hisat2.sam
```

```
@HD VN:1.0 SO:unsorted
@SQ SN:chr1 LN:195471971
@SQ SN:chr10 LN:130694993
@SQ SN:chr11 LN:122082543
@SQ SN:chr12 LN:120129022
@SQ SN:chr13 LN:120421639
@SQ SN:chr14 LN:124902244
@SQ SN:chr15 LN:104043685
@SQ SN:chr16 LN:98207768
@SQ SN:chr17 LN:94987271
```

Cabecera

```
(05MBIF) [UNIVERSIDADVIU\paula.soler@a-3edhijmgygwxr 03.Alineamiento]$ grep -v "^@" SRR1552444_hisat2.sam | head
```

[illegible]

Alineamiento

Archivo de alineamientos SAM (11 + 1)

SOLEXA-1GA-2_2_FC20EMB:5:251:979:328	0	chr1	10145	25	36M	*	0	0	AACCCCTAACCCCTAACCCCTAACCCCTAACT	hhhhHcWwhHTghcKA_ONhAAEEBZE?H?CBC?DA	NM:i:1	X1:
SOLEXA-1GA-2_2_FC20EMB:5:102:214:278	0	chr1	10148	25	36M	*	0	0	CCCTAACCCCTAACCCCTAACCCCTAACCTAAC	hbfbhhhXUYhT_ULZdLRTKNIMIKGLJCHFFJQJN	NM:i:0	X0:
SOLEXA-1GA-2_2_FC20EMB:5:195:284:685	16	chr1	10149	25	36M	*	0	0	CCAAACACTAACCCCTAACCCCTAACCCCTAAC	><>B@>?>?D>>?B?D>DBC?E@BDHAKCEKERL00	NM:i:1	X1:

QNAME: Query template name

FLAG: Combination of bit-wise flag

RNAME: Reference sequence name

POS: 1-based leftmost mapping position

MAPQ: Mapping Quality

CIGAR: CIGAR STRING

RNEXT: Reference name of the mate/next read

PNEXT: Position of the mate/next read

TLEN: Observed template length

SEQ: Segment sequence

QUAL: ASCII of Phred-scaled base quality

Optional metadata

Archivo de alineamientos SAM

```
SOLEXA-1GA-2_2_FC20EMB:5:251:979:328 0 chr1 10145 25 36M * 0 0 AACCCCTAACCTAACCTAACCTAACCTAACT hhhhHcWwhHTghcKA_ONhAAEEBZE?H?CBC?DA NM:i:1 X1:1
SOLEXA-1GA-2_2_FC20EMB:5:102:214:278 0 chr1 10148 25 36M * 0 0 CCCTAACCTAACCTAACCTAACCTAACCTAAC hbfbhhXUYhT_ULZdLRTKNIMIKGLJCHFFJQJN NM:i:0 X0:0
SOLEXA-1GA-2_2_FC20EMB:5:195:284:685 16 chr1 10149 25 36M * 0 0 CCAAACTAACCTAACCTAACCTAACCTAACCC ><>B@>?>?D>>?B?D>DBC?E@BDHAKCEKERL00 NM:i:1 X1:1
```

Optional metadata

QUAL:ASCII of Phred-scaled base quality

Bit	Description
1	0x1 template having multiple segments in sequencing
2	0x2 each segment properly aligned according to the aligner
4	0x4 segment unmapped
8	0x8 next segment in the template unmapped
16	0x10 SEQ being reverse complemented
32	0x20 SEQ of the next segment in the template being reverse complemented
64	0x40 the first segment in the template
128	0x80 the last segment in the template
256	0x100 secondary alignment
512	0x200 not passing filters, such as platform/vendor quality controls
1024	0x400 PCR or optical duplicate
2048	0x800 supplementary alignment


FLAG: Combination of bit-wise flag

QNAME: Query template name

Archivo de alineamientos SAM

```
SOLEXA-1GA-2_2_FC20EMB:5:251:979:328 0 chr1 10145 25 36M * 0 0 AACCCCTAACCTAACCTAACCTAACCTAACT hhhhHcWwhHTghcKA_ONhAAEEBZE?H?CBC?DA NM:i:1 X1:  
SOLEXA-1GA-2_2_FC20EMB:5:102:214:278 0 chr1 10148 25 36M * 0 0 CCCTAACCTAACCTAACCTAACCTAACCTAAC hbfbhhXUYhT_ULZdLRTKNIMIKGLJCHFFJQJN NM:i:0 X0:  
SOLEXA-1GA-2_2_FC20EMB:5:195:284:685 16 chr1 10149 25 36M * 0 0 CCAAACTAACCTAACCTAACCTAACCTAACCTAAC ><>B@>?>?D>>?B?D>DBC?E@BDHAKCEKERL00 NM:i:1 X1:
```

Una flag de 0 estaría relacionada con que la lectura se mapea en la cadena +, mientras que la flag de 16 indica que se está mapeando en al reverse



Archivo de alineamientos SAM

Picard
build passing

Latest Jar Release | Source Code ZIP File | Source Code TAR Ball | View On GitHub

A set of command line tools (in Java) for manipulating high-throughput sequencing (HTS) data and formats such as SAM/BAM/CRAM and VCF.

Decoding SAM flags

This utility makes it easy to identify what are the properties of a read based on its SAM flag value, or conversely, to find what the SAM Flag value would be for a given combination of properties.

To decode a given SAM flag value, just enter the number in the field below. The encoded properties will be listed under Summary below, to the right.

SAM Flag: **Explain**

Switch to mate Toggle first in pair / second in pair


Find SAM flag by property:
To find out what the SAM flag value would be for a given combination of properties, tick the boxes for those that you'd like to include. The flag value will be shown in the SAM Flag field above.

☐ read paired
☐ read mapped in proper pair
☐ read unmapped
☐ mate unmapped
☐ read reverse strand
☐ mate reverse strand
☐ first in pair
☐ second in pair
☐ not primary alignment
☐ read fails platform/vendor quality checks
☐ read is PCR or optical duplicate
☐ supplementary alignment

Summary:

<https://broadinstitute.github.io/picard/explain-flags.html>

Archivo de alineamientos SAM



Calidad de mapeo (MAPQ)

¿Qué valores podemos diferenciar de MAPQ en nuestro archivo de **SRR1552444_hisat2.sam** ?

```
grep -v "^@" SRR1552444_hisat2.sam | cut -f5 | sort -u
```

0 -> No se ha podido realizar una alineación única para esa lectura.

1 -> Alineamiento único de baja calidad.

60 -> Alineamiento que mapea de forma única.

Archivo de alineamientos SAM

```
SOLEXA-1GA-2_2_FC20EMB:5:251:979:328 0 chr1 10145 25 36M * 0 0 AACCCCTAACCCCTAACCCCTAACCCCTAACT hhhhHcWhhHTghcKA_ONhAAEEBZE?H?CBC?DA NM:i:1 X1:  
SOLEXA-1GA-2_2_FC20EMB:5:102:214:278 0 chr1 10148 25 36M * 0 0 CCCTAACCCCTAACCCCTAACCCCTAACCCCTAAC hbfbhhhXUYhT_ULZdLRTKNIMIKGLJCHFFJQJN NM:i:0 X0:  
SOLEXA-1GA-2_2_FC20EMB:5:195:284:685 16 chr1 10149 25 36M * 0 0 CCAAACTAACCCCTAACCCCTAACCCCTAACCCCTAAC ><>B@>?>?D>>?B?D>DBC?E@BDHAKCEKERL00 NM:i:1 X1:
```

RefPos: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19
Reference: C C A T A C T G A A C T G A C T A A C
Read: ACTAGAATGGCT

Aligning these two:

RefPos: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19
Reference: C C A T A C T G A A C T G A C T A A C
Read: A C T A G A A T G G C T

With the alignment above, you get:

POS: 5
CIGAR: 3M1I3M1D5M

CIGAR: CIGAR STRING

MAPQ: Mapping Quality

POS: 1-based leftmost mapping position


RNAME: Reference sequence name

FLAG: Combination of bit-wise flag

QNAME: Query template name

Esta codificación significa: que en la posición 5 han hecho match (M) los 3 primeros nucleótidos (empieza a alinear en la posición 5). En la siguiente hay 1 inserción, (1I) porque aparece en nuestra lectura pero no en el genoma de referencia, luego otras 3M, y luego 1D, es decir, una deletión, porque está en el genoma de referencia pero no en nuestra lectura


Archivo de alineamientos SAM



Compresión del archivo SAM a BAM



```
conda install bioconda::samtools  
conda install -c bioconda samtools
```



```
samtools view -Sbh SRR1552444_hisat2.sam > SRR1552444_hisat2.bam
```

```
samtools sort SRR1552444_hisat2.bam -o SRR1552444_hisat2.sorted.bam
```

```
samtools index SRR1552444_hisat2.sorted.bam
```

```
-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler users 2.7G Feb  6 11:14 SRR1552444_hisat2.bam  
-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler users 8.3G Feb  5 20:31 SRR155244_hisat2.sam  
-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler users 1.7G Feb  6 11:50 SRR155244_hisat2.sorted.bam  
-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler users 2.3M Feb  6 13:59 SRR155244_hisat2.sorted.bam.bai
```

Samtools stats

```
$ samtools stats SRR1552444_hisat2.sorted.bam > SRR155244_hisat2.sorted.bam.stats
```

```
$ cat SRR155244.sorted.bam.stats
```

```
# This file was produced by samtools stats (1.3.1+htslib-1.3.1) and can be plotted using plot-bamstats
```

```
# This file contains statistics for all reads.
```

```
# The command line was: stats SRR155244.sorted.bam
```

```
# CHK, Checksum [2]Read Names [3]Sequences [4]Qualities
```

```
# CHK, CRC32 of reads which passed filtering followed by addition (32bit overflow)
```

```
CHK 229871e0 5d4b8724 e32ebaf0
```

```
# Summary Numbers. Use `grep ^SN | cut -f 2-` to extract this part.
```

```
SN raw total sequences: 27906762
```

```
SN filtered sequences: 0
```

```
SN sequences: 27906762
```

```
SN is sorted: 1
```

```
SN 1st fragments: 27906762
```

```
SN last fragments: 0
```

```
SN reads mapped: 27089111
```

```
SN reads mapped and paired: 0 # paired-end technology bit set + both mates mapped
```

```
SN reads unmapped: 817651
```

```
SN reads properly paired: 0 # proper-pair bit set
```

```
SN reads paired: 0 # paired-end technology bit set
```

```
SN reads duplicated: 0 # PCR or optical duplicate bit set
```

```
SN reads MQ0: 89883 # mapped and MQ=0
```

```
SN reads QC failed: 0
```

```
SN non-primary alignments: 0
```



conda install
bioconda::rseqc



Table of Contents

RSeQC: An RNA-seq Quality Control Package

Release history

Download

- Download RSeQC
- Download test datasets
- Download gene models (update on 12/14/2021)
- Download ribosome RNA (update on 07/08/2015)

Installation

- Use pip3 to install RSeQC (v3.0.0 or newer)
- Use pip2 to install RSeQC (v2.6.6 or older)

Input format

Fetch chromosome size file from UCSC

Usage Information

- bam2fq.py
- bam2wig.py
- bam_stat.py
- clipping_profile.py
- deletion_profile.py
- divide_bam.py
- FPKM_count.py

RSeQC: An RNA-seq Quality Control Package

RSeQC package provides a number of useful modules that can comprehensively evaluate high throughput sequence data especially RNA-seq data. Some basic modules quickly inspect sequence quality, nucleotide composition bias, PCR bias and GC bias, while RNA-seq specific modules evaluate sequencing saturation, mapped reads distribution, coverage uniformity, strand specificity, transcript level RNA integrity etc.

Release history

RSeQC v5.0.1

- Oct 20, 2022
- Fix a bug in `scbam.py` to make it compatible with the latest `pysam` (v0.19.1).

RSeQC v5.0.0

- Oct 16, 2022
- add these functions to QC scRNA-seq data. * `sc_bamStat.py` * `sc_editMatrix.py` * `sc_seqLogo.py` * `sc_seqQual.py`

RSeQC v4.0.0

- Aug. 21, 2020
- Add `FPKM-UQ.py` to calculate HTSeq count, FPKM and `FPKM-UQ` values defined by TCGA
- `FPKM-UQ.py` could exactly reproduce TCGA FPKM-UQ values, if you use TCGA BAM file (or follow TCGA RNA-seq alignment workflow to generate your own BAM file), the `GDC.h38 GENCODE v22 GTF` file and the `GDC.h38 GENCODE TSV` file.

RSeQC v3.0.1

<https://rseqc.sourceforge.net/>

QualiMap



conda install
bioconda::qualimap

QualiMap

Evaluating next generation
sequencing alignment data

What is it?

Qualimap 2 is a platform-independent application written in Java and R that provides both a Graphical User Interface (GUI) and a command-line interface to facilitate the quality control of alignment sequencing data and its derivatives like feature counts.

Supported types of experiments include:

- Whole-genome sequencing
- Whole-exome sequencing
- RNA-seq (special mode available)
- ChIP-seq

How does it work?

Qualimap examines sequencing alignment data in **SAM/BAM** files according to the features of the mapped reads and provides an overall view of the data that helps to detect biases in the sequencing and/or mapping of the data and eases decision-making for further analysis.

Starting from version 2.0 QualiMap provides **multi-sample** comparison of alignment and counts data.

Features

- Fast analysis across the reference of genome coverage and nucleotide distribution;
- Easy to interpret summary of the main properties of the alignment data;
- Analysis of the reads mapped inside/outside of the regions provided in GFF format;
- Computation and analysis of read counts obtained from intersection of read alignments with genomic features;
- Analysis of the adequacy of the sequencing depth in RNA-seq experiments;
- Multi-sample comparison of alignment and counts data;
- Clustering of epigenomic profiles.

Download

Latest package for GNU Linux, MacOS and MS Windows:
[qualimap_v2.3.zip](#)

[Latest development snapshots](#)

[Public code repository](#)

[CIFP BioInfo local Maven repository](#)

[Version history](#)

Documentation

[QualiMap Online User Manual](#)

Sample data and output can be found [here](#)

Bioinformatics links:

Picard: a Java API (SAM-JDK) for creating programs that read and write SAM files.

FastQC: a quality control tool for high throughput sequence data.

SAMTools: essential utilities for manipulating alignments in the SAM format.

NOISeg: quality control and differential gene expression analysis for RNA-seq data.

Repitools: quality assessment, visualization, summarization and statistical analysis of epigenomics experiments.

Support


To get quick updates or report bugs/ suggestions please join the QualiMap Google group:



Your email:

<http://qualimap.conesalab.org/>

Flujo de trabajo del análisis de datos de RNA-seq (NGS)



GATK: Marcaje de duplicados




conda install
bioconda::picard

**Permite marcar y
eliminar duplicados,
aunque
normalmente se van
a mantener esos
duplicados.**

The screenshot shows the GATK documentation page for the MarkDuplicates tool. At the top, there's a navigation bar with links: User Guide, Tool Index, Blog, Forum, DRAGEN-GATK, Events, Download GATK4, and Sign in. Below this is a large blue banner with the text "Need Help?" and "Search our documentation", followed by a search input field. The main content area has a breadcrumb trail: GATK / Tool Index / 4.0.1.1. The title "MarkDuplicates (Picard)" is followed by a "Follow" button. Below the title is a profile icon for the "GATK Team" and a timestamp "September 19, 2022 11:23 · Updated". The description states: "Identifies duplicate reads. This tool locates and tags duplicate reads in a BAM or SAM file, where duplicate reads are defined as originating from a single fragment of DNA. Duplicates can arise during sample preparation e.g. library construction using PCR. See also [EstimateLibraryComplexity](#) for additional notes on PCR duplication artifacts. Duplicate reads can also result from a single amplification cluster, incorrectly detected as multiple clusters by the optical sensor of the sequencing instrument. These duplication artifacts are referred to as optical duplicates." On the right side, under the heading "Articles in this section", there is a list of related tools: "* Tool Documentation Index", "AccumulateVariantCallingMetrics (Picard)", "AddCommentsToBam (Picard)", "AddOATag (Picard)", and "AddOrReplaceReadGroups (Picard)".

<https://gatk.broadinstitute.org/hc/en-us/articles/360037052812-MarkDuplicates-Picard->

Flujo de trabajo del análisis de datos de RNA-seq (NGS)



Archivo de anotación de *Mus musculus* (formato GTF)



[Human](#) [Mouse](#) [How to access data](#) [FAQ](#) [Documentation](#) [About us](#)



Release M10 (GRCm38.p4)

- [Statistics of this release](#)
- [More information about this assembly](#) (including patches, scaffolds and haplotypes)

GTF / GFF3 files

Content	Regions	Description	Download
Comprehensive gene annotation	CHR	<ul style="list-style-type: none">• It contains the comprehensive gene annotation on the reference chromosomes only• This is the main annotation file for most users	GTF GFF3
Comprehensive gene annotation	ALL	<ul style="list-style-type: none">• It contains the comprehensive gene annotation on the reference chromosomes, scaffolds, assembly patches and alternate loci (haplotypes)• This is a superset of the main annotation file	GTF GFF3
Comprehensive gene annotation	PRI	<ul style="list-style-type: none">• It contains the comprehensive gene annotation on the primary assembly (chromosomes and scaffolds) sequence regions• This is a superset of the main annotation file	GTF GFF3
Basic gene annotation	CHR	<ul style="list-style-type: none">• It contains the basic gene annotation on the reference chromosomes only• This is a subset of the corresponding comprehensive annotation, including only those transcripts tagged as 'basic' in every gene	GTF GFF3
Basic gene annotation	ALL	<ul style="list-style-type: none">• It contains the basic gene annotation on the reference chromosomes, scaffolds, assembly patches and alternate loci (haplotypes)• This is a subset of the corresponding comprehensive annotation, including only those transcripts tagged as 'basic' in every gene	GTF GFF3
Long non-coding RNA gene annotation	CHR	<ul style="list-style-type: none">• It contains the comprehensive gene annotation of lncRNA genes on the reference chromosomes• This is a subset of the main annotation file	GTF GFF3
PolyA feature annotation	CHR	<ul style="list-style-type: none">• It contains the polyA features (polyA_signal, polyA_site, pseudo_polyA) manually annotated by HAVANA on the reference chromosomes• This dataset does not form part of the main annotation file	GTF GFF3

https://www.gencodegenes.org/mouse/release_M10.html

Archivo de anotación de *Mus musculus* (formato GTF)

```
(05MBIF) Data]$ tree -L 4
├── Annotation
├── Processed
│   ├── 01.Quality_control
│   │   ├── SRR1552444_fastqc.html
│   │   └── SRR1552444_fastqc.zip
│   ├── 02.Trimming
│   │   ├── SRR1552444.fastq.gz_trimming_report.txt
│   │   └── SRR1552444_trimmed.fq.gz
│   └── 03.Alignment
│       └── SRR1662444_hisat2.sam
├── Raw
│   └── SRR1552444.fastq.gz
└── Reference_genome
    ├── mm10
    │   ├── genome.1.ht2
    │   ├── genome.2.ht2
    │   ├── genome.3.ht2
    │   ├── genome.4.ht2
    │   ├── genome.5.ht2
    │   ├── genome.6.ht2
    │   ├── genome.7.ht2
    │   ├── genome.8.ht2
    │   ├── make_mm10.sh
    │   └── mm10_genome.tar.gz
```

```
$ gunzip gencode.vM10.annotation.gtf.gz
```

```
-rw-r--r-- 1 UNIVERSIDADVIU\paula.soler 765M May 1 2015
gencode.vM10.annotation.gtf.gz
```



General Feature Format (GFF3) / Gene Transfer format (GTF)

- Ambos formatos presentan en general la misma estructura.
- Ambos presentan 9 columnas de datos.
- Los campos se separan por tabulaciones (todos los campos deben contener un valor).

Sample GTF output from Ensembl data dump:

```
1 transcribed_unprocessed_pseudogene gene 11869 14409 . + . gene_id "ENSG00000223972"; gene_name "DDX11L1"; gene_source "havana"; gene_biotype "transcribed_unprocessed_pseudogene";
1 processed_transcript transcript 11869 14409 . + . gene_id "ENSG00000223972"; transcript_id "ENST00000456328"; gene_name "DDX11L1"; gene_source "havana"; gene_biotype "transcribed_unprocessed_transcript";
```

Sample GFF output from Ensembl export:


```
X Ensembl Repeat 2419108 2419128 42 . . hid=trf; hstart=1; hend=21
X Ensembl Repeat 2419108 2419410 2502 - . hid=AluSx; hstart=1; hend=303
X Ensembl Repeat 2419108 2419128 0 . . hid=dust; hstart=2419108; hend=2419128
X Ensembl Pred.trans. 2416676 2418760 450.19 - 2 genscan=GENSCAN00000019335
X Ensembl Variation 2413425 2413425 . + .
X Ensembl Variation 2413805 2413805 . + .
```

General Feature Format (GFF3) / Gene Transfer format (GTF)


```
(base) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr Annotation]$ ls
gencode.vM10.annotation.gtf
(base) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr Annotation]$ head gencode.vM10.annotation.gtf
##description: evidence-based annotation of the mouse genome (GRCm38), version M10 (Ensembl 85)
##provider: GENCODE
##contact: gencode-help@sanger.ac.uk
##format: gtf
##date: 2016-07-19
chr1 HAVANA gene 3073253 3074322 . + . gene_id "ENSMUSG000000102693.1"; gene_type "TEC"; gene_status "KNOWN"; gene_name "4933401J01Rik"; level 2; havana_gene "OTTMUSG000000049935.1";
chr1 HAVANA transcript 3073253 3074322 . + . gene_id "ENSMUSG000000102693.1"; transcript_id "ENSMUST000000193812.1"; gene_type "TEC"; gene_status "KNOWN"; gene_name "4933401J01Rik"; transcript_type "TEC"; transcri
pt_status "KNOWN"; transcript_name "4933401J01Rik-001"; level 2; transcript_support_level "NA"; tag "basic"; havana_gene "OTTMUSG000000049935.1"; havana_transcript "OTTMUST000000127109.1";
chr1 HAVANA exon 3073253 3074322 . + . gene_id "ENSMUSG000000102693.1"; transcript_id "ENSMUST000000193812.1"; gene_type "TEC"; gene_status "KNOWN"; gene_name "4933401J01Rik"; transcript_type "TEC"; transcript_statu
s "KNOWN"; transcript_name "4933401J01Rik-001"; exon_number 1; exon_id "ENSMUSE000001343744.1"; level 2; transcript_support_level "NA"; tag "basic"; havana_gene "OTTMUSG000000049935.1"; havana_transcript "OTTMUST000000127109.1";
chr1 ENSEMBL gene 3102016 3102125 . + . gene_id "ENSMUSG000000064842.1"; gene_type "snRNA"; gene_status "KNOWN"; gene_name "Gm26206"; level 3;
chr1 ENSEMBL transcript 3102016 3102125 . + . gene_id "ENSMUSG000000064842.1"; transcript_id "ENSMUST000000082908.1"; gene_type "snRNA"; gene_status "KNOWN"; gene_name "Gm26206"; transcript_type "snRNA"; transcript
status "KNOWN"; transcript_name "Gm26206-201"; level 3; transcript_support_level "NA"; tag "basic";
```

Columna	Tipo	Descripción
1	SEQNAME	Nombre del cromosoma o <i>scaffold</i> .
2	SOURCE	Nombre del programa de predicción
3	FEATURE	Categoría de secuencia: <i>CDS</i> , <i>intron</i> , <i>exon</i> , <i>gene</i> , etc.
4	START	Inicio
5	END	Fin
6	SCORE	Vacío normalmente “.”
7	STRAND	Cadena (+) o (-) dónde se encuentra la .
8	FRAME	Vacío normalmente “.”
9	ATTRIBUTE	Información adicional separada por “;”.

Flujo de trabajo del análisis de datos de RNA-seq (NGS)




Recuento del número de lecturas por gen anotado




Exon	Number of reads
gene1 - exon1	3
gene1 - exon2	2
gene2 - exon1	3
gene2 - exon2	4
gene2 - exon3	3

Recuento del número de lecturas por gen anotado



- Número de lecturas por exón
- Número final de lecturas por gen, entiendo que este es la unión de todos los exones.

Recuento del número de lecturas por gen anotado



HTSeq-count (Anders *et al.* 2015)

featureCounts (Liao *et al.* 2013)

Recuento de los alineamientos (HTseq-count)

HTSeq 0.11.1 documentation »

Table Of Contents

Counting reads in features with `htseq-count`

- Usage
- Options
- Frequently asked questions

Previous topic

Quality Assessment with `htseq-ga`

Next topic

Version history

This Page

Show Source

Quick search

Go

Counting reads in features with `htseq-count`

Given a file with aligned sequencing reads and a list of genomic features, a common task is to count how many reads map to each feature.

A feature is here an interval (i.e., a range of positions) on a chromosome or a union of such intervals.

In the case of RNA-Seq, the features are typically genes, where each gene is considered here as the union of all its exons. One may also consider each exon as a feature, e.g., in order to check for alternative splicing. For comparative ChIP-Seq, the features might be binding region from a pre-determined list.

Special care must be taken to decide how to deal with reads that align to or overlap with more than one feature. The `htseq-count` script allows to choose between three modes. Of course, if none of these fits your needs, you can write your own script with HTSeq. See the chapter [A tour through HTSeq](#) for a step-by-step guide on how to do so. See also the FAQ at the end, if the following explanation seems too technical.

The three overlap resolution modes of `htseq-count` work as follows. For each position i in the read, a set $S(i)$ is defined as the set of all features overlapping position i . Then, consider the set S , which is (with i running through all position within the read or a read pair)

- the union of all the sets $S(i)$ for mode `union`. This mode is recommended for most use cases.
- the intersection of all the sets $S(i)$ for mode `intersection-strict`.
- the intersection of all non-empty sets $S(i)$ for mode `intersection-nonempty`.

If S contains precisely one feature, the read (or read pair) is counted for this feature. If S is empty, the read (or read pair) is counted as `no_feature`. If S contains more than one feature, `htseq-count` behaves differently based on the `--nonunique` option:

- `--nonunique none` (default): the read (or read pair) is counted as `ambiguous` and not counted for any features. Also, if the read (or read pair) aligns to more than one location in the reference, it is scored as `alignment_not_unique`.
- `--nonunique all`: the read (or read pair) is counted as `ambiguous` and is also counted in all features to which it was assigned. Also, if the read (or read pair) aligns to more than one location in the reference, it is scored as `alignment_not_unique` and also separately for each location.

Notice that when using `--nonunique all` the sum of all counts will not be equal to the number of reads (or read pairs), because those with multiple alignments or overlaps get scored multiple times.

The following figure illustrates the effect of these three modes and the `--nonunique` option:

https://htseq.readthedocs.io/en/release_0.11.1/count.html

Recuento de los alineamientos (HTseq-count).

¿Qué contar?

⇒ *id_attr=gene_id*

Para recuentos por gen



¿Dónde contar?







⇒ *tipo=exón*

Contar alineaciones superpuestas exones

¿Como contar?

⇒ *modo=unión*

Para la expresión diferencial en el nivel genético

	union	intersection_strict	intersection_nonempty
	gene_A	gene_A	gene_A
	gene_A	no_feature	gene_A
	gene_A	no_feature	gene_A
	gene_A	gene_A	gene_A
	gene_A	gene_A	gene_A
	ambiguous	gene_A	gene_A
	ambiguous	ambiguous	ambiguous

Recuento de los alineamientos (HTseq-count).



```
conda install bioconda::htseq
```

```
(05MBIF) $ htseq- (+TAB)
```

```
htseq-count      htseq-count-barcodes  htseq-qa
```

con -t le decimos que busque los exones, -i que agrupe las lecturas que pertenecen a un mismo gen con el id, --stranded le estamos diciendo que NO tiene información de hebra específica, -f el formato del archivo, -r que lo ordene por la posición, -s va junto con la opción stranded, que básicamente es para especificar que no tenemos información de orientación específica.

<https://anaconda.org/bioconda/htseq>

Lanzar la instrucción

```
(05MBIF) [03.Alignment]$ htseq-count -t exon -i gene_id --stranded=no -f bam -r pos -s no SRR1552444_hisat2.sorted.bam  
../Annotation/gencode.vM10.annotation.gtf > ../Results/SRR1552444_counts.tsv
```

100000 GFF lines processed.

200000 GFF lines processed.

300000 GFF lines processed.

Archivo de recuentos (formato tsv)

```
(05MBIF) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr Results]$ head SRR1552444_counts.tsv
```

```
ENSMUSG000000000001      5440
```

```
ENSMUSG000000000003      0
```

```
ENSMUSG000000000028     256
```

```
ENSMUSG000000000031     109
```

```
ENSMUSG000000000037      55
```

```
ENSMUSG000000000049      1
```

```
ENSMUSG000000000056     266
```

```
ENSMUSG000000000058     2800
```

```
ENSMUSG000000000078    45371
```

La primera columna es un identificador de los genes recogidos en este archivo y la segunda columna es el número de lecturas asociadas a ese gen

Archivo de recuentos (formato tsv)

```
(base) [UNIVERSIDADVIU\paula.soler@a-3edhijmqygwxr Results]$ tail SRR1552444_counts.tsv
```


ENSMUSG00000106667	0
ENSMUSG00000106668	0
ENSMUSG00000106669	0
ENSMUSG00000106670	0
ENSMUSG00000106671	1

__no_feature	2342826	número de lecturas que no se han podido asociar a ninguna característica genómica
__ambiguous	1002121	lecturas que se han podido alinear a una o más de una característica genómica, una lectura que se ha superpuesto en dos genes y no se a qué gen pertenece y por tanto genera ambigüedad
__too_low_aQual	3261453	lecturas que no se contaron porque la calidad de su alineamiento era bajo
__not_aligned	817651	numero de lecturas que no se pudieron alinear contra el genoma de referencia
__alignment_not_unique	0	número de lecturas que tienen más de un alineamiento

Merging de todas las muestras obtenidas: *Count matrix*

samples: want to see if differences across condition are significant
(w.r.t. biological and technical variation)

features (e.g. genes)



	SRR1039508	SRR1039509	SRR1039512	SRR1039513	SRR1039516
ENSG000000000003	679	448	873	408	1138
ENSG000000000005	0	0	0	0	0
ENSG000000000419	467	515	621	365	587
ENSG000000000457	260	211	263	164	245
ENSG000000000460	60	55	40	35	78

Package "EdgeR"



viu

Universidad
Internacional
de Valencia